CAMERA

SWIFTCAM - DIGITAL VIDEO CAMERA 1-CCD

- Multilingual on-screen-menu
- Digital zoom
- Camera control with keyboard
- Character generator



CAMERA

SWIFTCAM - DIGITAL VIDEO CAMERA 1-CCD

Technical Data

Camera head: • Aspect ratio: 4:3

· Image sensor: 1/3" CCD

· Immersible (using the attached cap)

Voltage range (current consumption): • 100-240V (25 VA)

Frequency: • 50-60 Hz
Device protection: • 2x T3,15 A

Video system: • PAL/NTSC

 \cdot Image resolution: PAL 752 x 582/ NTSC 768 x 494

TV resolution: > 470 linesSensitivity: < 1 lux (1/1.4)

· White balance automatically, with data hold

shutter speed PAL (NTSC): 1/50 (1/60) to 1/100.000 sec

Video out terminal: 2x Video; 2x S-Video; 1x Y-Pb-Pr (Component); 1x SDI; 1x USB 2.0

Protection class (acc.to IEC 60601-1): I Safety class: BF IP code: 20

Standard: • IEC 60601-1-2

Device classification (acc. to MDD 93/42/EEC): \cdot I

Dimensions (WxHxD): 355mm x 103mm x 253mm

Weight: • 4,4 kg

Delivery contains • CCU, camera head with 3,05m (10ft) cable, mains cable, User manual,

Y/C cable, BNC cable, HDMI - DVI cable (optional), HDMI cable (optional)

Series	SwiftCAM	SwiftCAM-D	
On-screen-display	X	Х	
Menu	Х	Х	
Titel generator option PS/2	X	X	
HDMI out	-	X	
Zoom function	Х	Х	

Options

Freeze function: • Freeze the image via push buttons

Remote camera head: • With programmable push buttons (zoom, freeze)

Remote footswitch: • With programmable push buttons

SwiftCAM camera head



Weight: 80 g Size: d= 29 mm; l= 55 mm



Weight: 105 g

Size: d= 34 mm; l= 73 mm

Without push buttons