


Prismxel

Technical Specification

Model	Prismxel						
Laser Type	CO2 Laser						
Wave length	10.600nm						
Laser power	1- 40W						
Guide Beam	650nm Less than 5mW						
Function							
Aiming Beam	Diode Laser 655nm (5mw) :1 - 10 Step						
Handpiece Type	50mm,100mm(Basic) E.N.T & Dental Handpiece Optional						
Surgical Mode	<table border="1"><tr><td>Laser Operation Mode</td><td>1.CW 2.Normal Pulse 3.Super Pulse 4.Super Pulse Continuous 5.Ultra Pulse 6.Smart Protocol(General Surgery/Dental Surgery)</td></tr><tr><td>Pulse Duration</td><td>• Ultra Pulse :90us- 900us • Super Pulse :1000us - 3500us • Normal Pulse :4ms - 20ms</td></tr><tr><td>Reapet Time</td><td>OFF,2ms,5ms,10ms,30ms,50ms,70ms,100ms,250ms,500ms</td></tr></table>	Laser Operation Mode	1.CW 2.Normal Pulse 3.Super Pulse 4.Super Pulse Continuous 5.Ultra Pulse 6.Smart Protocol(General Surgery/Dental Surgery)	Pulse Duration	• Ultra Pulse :90us- 900us • Super Pulse :1000us - 3500us • Normal Pulse :4ms - 20ms	Reapet Time	OFF,2ms,5ms,10ms,30ms,50ms,70ms,100ms,250ms,500ms
Laser Operation Mode	1.CW 2.Normal Pulse 3.Super Pulse 4.Super Pulse Continuous 5.Ultra Pulse 6.Smart Protocol(General Surgery/Dental Surgery)						
Pulse Duration	• Ultra Pulse :90us- 900us • Super Pulse :1000us - 3500us • Normal Pulse :4ms - 20ms						
Reapet Time	OFF,2ms,5ms,10ms,30ms,50ms,70ms,100ms,250ms,500ms						
General Features							
Electrical Requirements	220v, 6A. single phase,50/60hz						
Cooling System	water cooling						
Dimensions	350x360x1200						
Weight	45Kg						
Display	10.2 inches colorcd						

CE
0068



Prismxel



Prismxel → Surgical Mode/Dental Mode

Smart mode for Surgical Field (General Surgery , Dental Surgery)



General Surgery Indication

Corns, Warts etc	Freckles
Birthmarks(mole)	Skin resurfacing
Phimosectomy	Tumor elimination
Hemorrhoids and anal fistulas	Fine wrinkle
Condyloma	Syringoma
Snoring	Allergic splenitis

Dental Surgery Indication

The second implant operation	Hyper sensitive
Hemostasis after dental extraction	Herpetic gingivitis
Crown lengthening	Frenulum excision
Gingival excision, Gingival plastic operation	Cleaning of dental cuticle cyst
Removal of pigmentation	Pulp capping -Fluorine coating
Sterilization in root canal	Sterilization in parodontal cyst
Abscess treatment	

→ General Mode

General CO2 Laser Mode



General Surgery Indication

Corns, Warts etc
Birthmarks(mole)
Phimosectomy
Hemorrhoids and anal fistulas
Condyloma
Snoring
Freckles
Skin resurfacing
Tumor elimination
Fine wrinkle
Syringoma
Allergic splenitis