


RP11A

FLEXIBLE VIDEO RHINOLARYNGOSCOPE



160g

Ultra-light handle, more stable grip, reduce the fatigue of the doctor, and make the examination smoother.



1280x720P

3.2mm high-definition video endoscope, ultra-thin insertion, full of details

A variety of connection methods are suitable for various diagnosis and treatment environments.



Computer



Pad



Processor

Optical system

- | | |
|---------------------|---------|
| ■ Field of view | 100° |
| ■ Depth of field | 3-150mm |
| ■ Direction of view | 0° |

Insertion portion

- | | |
|---------------------|-----------------|
| ■ Insertion tube OD | 3.2mm/3.8mm |
| ■ Angulation range | U 130° / D 130° |
| ■ Working length | 330mm±10% |
| ■ Extension cable | USB、REMO |
| ■ Max. resolution | 1280x720P |



weseemed
medical technology

ZHONGSHAN WESEE MEDITECH CO., LTD.

6F, Block C, 2Jianshuo Str. HoChung Village,
Nanlang, Zhongshan 528451 P. R. China

TEL:(+86)0760-86960005

Email:info@weseemed.com

www.weseemed.com

