Lighted Articulating Ear Curette[™]

Bionix[®] Lighted Articulating Curette combines illumination. magnification and a 90° articulating tip to help healthcare providers take less time to see and safely remove ear canal passage obstructions.



PROBLEM

Cerumen plugs and foreign bodies in the ear are often difficult to remove. Practitioners lack the ability to visualize or work from behind the occlusion to remove it.

SOLUTION

The Lighted Articulating Ear Curette brilliantly illuminates the ear canal and provides an articulating tip that allows the practitioner to get behind the occlusion.

RESULTS

Direct visualization and an articulating ear curette eliminates the blind, top down sweeping motion that can abrade the ear canal improving safety and efficiency.

Refer to the package insert supplied with product for specific information on indications for use, contraindications, warnings, and precautions,

POTENTIAL REIMBURSEMENT (U.S. ONLY)

CPT Code/Description 69210 - For potential reimbursement of impacted cerumen. 69200 - For potential reimbursement of foreign body removal from the ear. 30300 - For potential reimbursement of foreign body removal from the nose.







Boulevard Général Wahis 53 • 1030 Brussels, BELGIUM Tel: +(32) 2.732.59.54 • Fax: +(32) 2.732.60.03 • E-Mail: mail@obelis.net Bionix LLC • 1670 Indian Wood Circle, Maumee, Ohio 43537, United States Phone: +1 419.727.8421 • Fax: +1 800.455.5678 • www.Bionix.com



PRODUCT OFFERING

- #2511 Lighted Articulating VersaScoop® (3mm tip) Includes 25 Articulating Ear Curette Tips, 1 Light Source and 1 Magnification Lens.
- #2200 Extra Light Source #2265 - Pack of 5 Magnification Lenses



INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS

800.551.7096 · www.bionix.com