# Lighted Forceps for Foreign Body Removal

Bionix® Lighted Forceps combine illumination and magnification to help healthcare providers take less time to see and safely remove ear canal and nasal passage obstructions.





## **PROBLEM**

Objects can be difficult to grasp or remove from the ear canal or nasal cavity when they are soft and organic in shape, such as leafy objects.

## **SOLUTION**

The Lighted Forceps provides illumination, magnification and incorporates the grasping ability of forceps for foreign body removal.

## **RESULTS**

An intuitive, easy to use tool that is great for removing foreign bodies of the ear and nose.

## PRODUCT OFFERING

#2750 - Box of 10 Forceps, 1 Light Source, 1 Magnification Lens

#2775 - Box of 10 Replacement

#2200 - Extra Light Source

#2265 - Pack of 5 Magnification Lenses

"I used the (Lighted) Forceps to remove a cotton wad from a woman's ear and it was sooooooo easy." - Judith Kuchta, ARNP

Refer to the package insert supplied with product for specific information on indications for use. contraindications, warnings, and precautions.

## POTENTIAL REIMBURSEMENT (U.S. ONLY)

CPT Code/Description

69200 - For potential reimbursement of foreign body removal from the ear. **30300** - For potential reimbursement of foreign body removal from the nose.















Bionix LLC • 1670 Indian Wood Circle, Maumee, Ohio 43537, United States Phone: +1 419.727.8421 • Fax: +1 800.455.5678 • www.Bionix.com



INNOVATIVE SOLUTIONS FOR EVERYDAY PROBLEMS